



Calf Puller For Difficult Calving

RMI-1002

Easy To Operate Zigzag Action, extra long 65" with 12" extension for a total size



Ratchet Style Calf Puller

RMI-1012

Easy to use zigzag action. Breech is supplied with an adjustable user-friendly ch



Puppy OB Forceps

RMI-1001

Stainless Steel
8"



Bovine Cervical Forceps

RMI-1003

Stainless Steel
12"



Krey Hook

RMI-1005

Stainless Steel Double Action OB Hook bites into part of Fetus in the process of embryotomy

3"



OB Wire Guide

RMI-1008

Stainless Steel It helps to bring wire around any tight or difficult location on the fetus

7"



Lindse's Embryotomy Knife

RMI-1004

Stainless Steel With one Large Finger tip and handle adapted to the palm
5"



Ostertag's Blunt Eye Hook

RMI-1006

Stainless Steel Hole Diameter 5/8" Hook Size 3/4"
2"



Vienna Double Blunt Eye Hook

RMI-1007

Stainless Steel Hinged for easier placement
5 1/2"



Suture Style OB Handle

RMI-1009

Knurled non-slip Grip
4"



Lyss OB Handle

RMI-1010

Stainless Steel Grasp to make it easy to change position on OB Wire
3"



Pig Obstetrical Forceps

RMI-1011

Stainless Steel Danish Model Forceps Open sideways. Softly hold piglet either aro
23"



Teat Tumor Extractor

RMI-1102

Stainless Steel with Two Cutting Edges



Teat Slitter

RMI-1101

Stainless Steel Blades are opened by loosening. Small Burr on Side of the Instrument
7 1/2"



Milk Tube Adjustable

RMI-1104

3"



Milking Tube

RMI-1103

With Olive Hub, Nickel Plated
2.0mm 2.5mm 3.0mm 4.0mm



Milk Tube

RMI-1105

Stainless Steel
2" 3" 1.5" 2.5"



Bayer Mouth Wedges

RMI-1301

A very useful tool to examine the oral cavity



Equine Speculum By Schoupe

RMI-1302

Stainless Steel Minor risk of slipping. Inserted without difficulty. Greater Free



Equine Molar Cutter

RMI-1305

Handle 11 1/2" Powerful combined action. Breakdown Model for Easy Storage



Equine Molar Forceps

RMI-1306

Stainless Steel



Thompson's Universal Dental Float

RMI-1303

Long Lasting Tungsten Carbide Side to side blade. 20" Shaft



Equivet Nose Twitch

RMI-1311

1m Nylon Rope and Snap Hook



Automatic Bull Lead

RMI-1207

Available in Brass,C/P and Aluminum



Harm's Bull Holder

RMI-1208

Stainless Steel



Spring Bull Lead

RMI-1209

These Spring lock leads are built for quick insertion and removal.They are place



Pig Holder With Lock

RMI-1210



Canine Mouth Gag

RMI-1211

Stainless Steel Spring Action with Nylon Inserts,available in three Sizes

4.5"5 3/4"6 3/4"



Hog Holder

RMI-1201

Small Especially Designed for handling small Swines



Cattle Twitch 17"

RMI-1202

Nylon Coated Handles for sure Grip



Killian Vaginal Speculum For Gogs

RMI-1203

Has Set Screw and 3 1/2"Blades



Polansky Vaginal Speculum With 3 Blades

RMI-1204

Stainless Steel Has Butterfly Screw with Two Blades on Top



Anti Kick Device

RMI-1205

Born Screw Down Model For Use on Achilles Tendon



Anti Kick Bar

RMI-1206

Excellent Quality Adjustable with Spring Lock Prevents Kicking at Milking and Under Treatment



Hoof Knife

RMI-1409

Stainless Steel High Quality Instrument with small
Curl Blade Width 5/8" Right or



Hoof Knife

RMI-1410

Stainless Steel Same as 1409 except featuring
Narrow 3/8" blade Right or Left Cu



Large Hoof Tester

RMI-1406

Stainless Steel Heavy Duty
18"



Hoof Test

RMI-1405

Stainless Steel
13"



Swiss Hoof Knife

RMI-1401

Stainless Steel Swiss Pattern Hard Wood Handle



Folding Hoof Knife

RMI-1402

Stainless steel Combined with Hoof Pick Handy Pocket Size



Hoof Pick

RMI-1403

Steel Hoof Pick With Vinyl Grip Allows you to Remove Foreign objects from your horses hoof



Foot Rot Hoof Shear

RMI-1404

Tempered Steel Available with Serrated and Non-serrated Blades Heavy Handles wit



Hoof Trimmer

RMI-1407

Strong Compound Action 20"Long



Hoof Knife

RMI-1408

Stainless Steel With Small Curl Double Edge Blade Width 5/8" Hard Wood Handle



Equivet Full Mouth Speculum

RMI-1312

Stainless Steel



Woolfthooth Extractor Forceps

RMI-1307

Stainless Steel

11 1/2"



REYNOLDS Tooth Extractor

RMI-1308

Stainless Steel Lower jaws

37 cm



House Combination Dental Float

RMI-1304

Stainless Steel Standard dental float with both straight and angled shaft with o

20"



REYNOLDS Tooth Extractor

RMI-1309

Stainless Steel Upper Jaws

37 cm



Baling Gun

RMI-1310

Nickel Plated Metal Diameter 30mmx

51 cm



Barne's Dehorner

RMI-1702

Metal Handles 13", 17"



Barne's Dehorner

RMI-1701

Wooden Handles 13"



Baling Gun

RMI-1703

With Spring Clips (1" diameter)



Baling Gun

RMI-1704

With Spring Clips (5/8" diameter)



Baling Gun

RMI-1705

Plastic Head With Spring Clip (1" diameter)



Baling Gun

RMI-1706

Plastic Hand With Spring Clip (5/8" diameter)



White's Modified Emasculator

RMI-1601

Stainless Steel Double Crush with extra crushing area
12"9"



Hausman Standard Emasculator

RMI-1603

Stainless Steel
12"7 1/2"9"



Reimer Emasculator

RMI-1604

Stainless Steel Double Crushing Part.The Cutting Lever is Worked Separately
12 1/2"



Serra Modified Emasculator

RMI-1605

Stainless Steel For One hand Use.Reduces Tissues Slippage and post-op hemorrhage
14"



Hetzel Emasculator

RMI-1606

Stainless Steel
5 1/2"



Laszlo Emasculator

RMI-1607

Stainless Steel
6 1/2"



Sand's Castrator
RMI-1602
 Stainless Steel Crushing Only



Castrating Plier
RMI-1608
 Stainless Steel



Castrating Plier
RMI-1609
 C/P



Castrating Plier
RMI-1610
 Lightweight, alternative to metal castrators Resistant and durable



Deluxe Castrating Plier
RMI-1611
 Aluminum



Emasculator For Bloodless Castration 9"
RMI-1612
 Forged Carbon Steel Nickle Plated



Rochester Pean

RMI-1807

Straight & Curved, Stainless Steel Also available with Tungsten Carbide
Size: 5.5", 6 1/4", 8"



Allis Tissue Forceps 6" & 7"

RMI-1808

Stainless Steel Also available with Tungsten Carbide



Kocher Forceps 1x2 Teeth

RMI-1809

Straight & Curved, Stainless Steel Also available with Tungsten Carbide
Size: 5.5", 6 1/4", 8"



Jones Artery Forceps 5"

RMI-1810

Stainless Steel Used in Dehorning in pulling the Corneal Artery



Jones Towel Clamps 3 1/2"

RMI-1811

Stainless Steel



Backhaus Towel Clamps 3 1/2"

RMI-1812

Stainless Steel



Halsted Mosquito Hemostats 5"

RMI-1801

Straight,Stainless Steel Also available with Tungsten Carbide



Halsted Mosquito Hemostats 5"

RMI-1802

Curved Stainless Steel Also available with tungsten Carbide



Kelly Hemostats 5 1/2"

RMI-1803

Straight,Stainless Steel Also available with Tungsten Carbide



Kelly Hemostats 5 1/2"

RMI-1804

Curved,Stainless Steel Also available with Tungsten Carbide



Circle Hemostats 6"

RMI-1805

Straight,Stainless Steel Also available with Tungsten Carbide



Circle Hemostats 6"

RMI-1806

Curved,Stainless Steel Also Available with Tungsten Carbide



Baling Gun

RMI-1707

Enlarged to give the foreign body magnets,flexible
Inner diameter 1.5"



Multi Bolus Baling Gun

RMI-1708



British Capsule Forceps

RMI-1709

Stainless Steel



One Hand Sheep Clipper

RMI-1710



Two Hand Sheep Clipper

RMI-1711



Black Blue Sheep Shear, 13.6"

RMI-1712

Made of Hardened and Tempered Steel.Rubber
Stoppers Smoothen Shearing. Straight



Flat S.S. Piglet

RMI-2004

FLAT, stainless steel, 1/2" thread
1/2"



Hex Brass

RMI-2005

HEX, Brass, 1/2" thread
1/2"



Hex Brass Piglet

RMI-2006

HEX, Brass, Piglet, 1/2" thread
1/2"



Hex S.S.

RMI-2007

HEX, stainless steel, 1/2" thread
1/2"



Hex S.S. Piglet

RMI-2008

HEX, stainless steel, Piglet, 1/2" thread
1/2"



Hex S.S. Piglet Narrow

RMI-2009

HEX, stainless steel, Piglet, Narrow, 1/2" thread
1/2"



Ear Tag Pliers

RMI-794

9.5"



Sheep Shears

RMI-1901

5 1/2" long

5 12"



Sheep Shears

RMI-1902

6 1/2" long, red handle

5 12"



Flat Brass

RMI-2001

FLAT, Brass, 1/2" thread

1/2"



Flat Brass Piglet

RMI-2002

FLAT, Brass, Piglet, 1/2" thread

5 12"



FLAT

RMI-2003

FLAT, stainless steel, 1/2" thread

1/2"



Meat Chopper

RMI-795



Meat Chopper

RMI-796



Sheep Shears

RMI-1905

Stainless Steel



Sheep Shears

RMI-1903

Stainless Steel



Sheep Shears

RMI-1904

Stainless Steel



Ear Tag Pliers

RMI-793

9.5"



Eldonain Peacock Safety Irons

RMI-3407

4"3 1/2"3 3/4"4 1/2"4 1/4"



Shires Bent Leg Stirrup Irons

RMI-3408

4"4 1/4" 4 3/4"3 3/4"4 1/2"



ProRMile Stirrups

RMI-3404



Stirrups Irons

RMI-3401

Stainless steel stirrup irons

5"4"4.5"3.5"4.75"4.25"



Korsteel Flexi Stainless Steel Stirrups

RMI-3402

4 1/4" 4 3/4"



Foot Saver Stirrups

RMI-3403

4"4 1/4" 4 3/4"



Premium Profile Stirrups

RMI-3405

AdultChilds



Premium Profile Stirrups

RMI-3406

AdultChilds





Hex S.S. Think Pin

RMI-2010

WET, HEX, stainless steel, think pin 1/2" thread
1/2"



Wet Hex S.S. Long Pin

RMI-2011

WET, HEX, stainless steel, long pin, 1/2" thread
1/2"



Pig Tooth Nipper With

RMI-1505

Stainless Steel

5 1/2"



Pig Tooth Nipper

RMI-1506

Stainless Steel With Roller Spring

5 1/2"



Veterinary Syringe

RMI-4301

Veterinary plastic Steel Syringe A-Type. 20ml. 1.
Material: poly 4-methyl-penten



Veterinary Syringe

RMI-4302

Veterinary plastic Steel Syringe B-Type. 20ml. 1.
Material: poly 4-methyl-penten



Veterinary Syringe

RMI-4303

Veterinary plastic Steel Syringe C-Type, 10ml, 20ml,
30ml, 50ml 1. Material: pol



Veterinary Syringe

RMI-4304

Veterinary plastic Steel Syringe D-Type, 10ml, 20ml,
30ml, 50ml 1. Material: pol



Veterinary Syringe

RMI-4305

Veterinary plastic Steel Syringe E-Type, 10ml, 20ml,
30ml, 50ml 1. Material: pol



Veterinary Syringe

RMI-4306

Veterinary plastic Steel Syringe F-Type, 10ml, 20ml,
30ml, 50ml 1. Material: pol



Ring-O-Matic Hog Catcher

RMI-4001

Made for one person use, the Ring Hog Catcher has a smooth, sanitary cable on th



Half Round Hoof Nipper

RMI-4101

The smallest of the nippers. Ideal for foals and any of the smaller-hoofed horses

12"



Hoof Nipper

RMI-4102

The 14" nipper is considered the workhorse of the nippers. It is especially usefu

14"





Snake Tong

RMI-2801

Gentle Giant Collapsible Snake Tong
25"40"52"60"72"



Snake Tong

RMI-2802

Gentle Snake Tong Gentle jaws made with fine
aluminum
25"40"52"60"72"



Snake Tong

RMI-2803

Snake Tongs Strong shaft to handle big reptiles
25"40"52"60"72"



Snake Tong

RMI-2804

Snake Tongs Strong shaft to handle big reptiles
25"40"52"60"72"



Snake Tong

RMI-2805

Snake Tongs Strong shaft to handle big reptiles
25"40"52"60"72"



Snake Tong

RMI-2806

Snake Tongs Strong shaft to handle big reptiles
25"40"52"60"72"



Either A Gelpi Or Weitlaner Retractors

[RMI-789](#)

7"



Bone Rongers Still

RMI-790

9"



Crinipel Wire Cutter

RMI-786

9"



Crinipel Wire Cutter

RMI-787

9"



Compound Action Crimper

RMI-788

7"



Crimping Forceps

RMI-791

8"



Stifle Retractors

RMI-792

8"



Barrel Dee Ring Snaffle

RMI-3202

These bits are made with 18/8 stainless steel The Dee Ring bit helps keep the bi
4"4.75"



Full Check Snaffle

RMI-3101

Full Cheeks aid in communication from the rider to pony by acting on the sides of the mouth.

4"4.5"3.5"



Full Check Eggbutt Snaffle

RMI-3102

These full cheek snaffles are made with 18/8 stainless steel. The cheek measurement is 4 inches.

4"4.75"



John Patterson Eggbutt Full Cheek

RMI-3103

The full cheek gives control on the sides of the mouth, the lips and the corners of the mouth.

4"4.5"4.75"



John Patterson Eggbutt Full Cheek

RMI-3104

The full cheek gives control on the sides of the mouth, the lips and the corners of the mouth.

4"4.5"4.75"



Snake Tong

RMI-2807

The jaw and handle of this tong is made from best quality aluminum with nice grip.

25"40"52"60"72"



Snake Hook

RMI-2901

Traditional snake hook stainless steel material

100 cm39.99"



Snake Hook

RMI-2902

Professional Snake Hook

45"



Snake Hook

RMI-2904

Mini Snake Hook For Small Snakes

24"



Snake Hook

RMI-2903

Standard Snake Hook Stainless Steel

40"



Snake Safety Gloves

RMI-3001

Snake Handling Gloves

SmallMediumLarge



Herm Sprenger Turnado Loose Ring Snaffle

RMI-3301

Herm Sprenger bits have outstanding German quality and are professionally developed.
4.25" With 70 mm (2.75") rings



Herm Sprenger KK Ultra Double Joined Bradoon

RMI-3302

Herm Sprenger bits have outstanding German quality and are professionally developed.
4.5" 4.125"



Polansky Veterinary Vaginal Speculum

RMI-3801

The speculum has 3 blades. Has butterfly screws with 2 Blades on top for exposure.

11"



Ring Applicators

RMI-3901

Castration Ring Applicator

9"



Ring Applicators

RMI-3903

Bull Ring Applicator

10"



Ring Applicators

RMI-3904

Castration Ring Applicator

9"